# Ella-Simple ELISA 系统

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| Description                    | Specification  |
|--------------------------------|--|
| Environmental Conditions       | For Indoor Use Only<br>Altitude up to 2000 m (6600 ft)<br>Temperature Range: 15° to 30° C (59° to 86° F)<br>Humidity: 15 to 80% RH, non-condensing<br>Rated Pollution Degree: Degree 2   |
| Electrical Ratings             | 100-240 V(AC), 300 VA, 45/65 Hz, Class I<br>Installation Category (Over Voltage) II  |
| Electrical Test Specifications | CB Scheme:<br>IEC 61010-1:2010<br>Supplemented by IEC 60825-1:2014<br><br>TÜV SÜD "CUE" Scheme:<br>CAN/CSA-C22.2 No. 61010-1:2012<br>UL 61010-1:2012<br>EN 61010-1:2010<br>EN 60825-1:2014   |
| Physical Conditions            | Overall dimensions:<br>37 cm (w) x 54 cm (d) x 26 cm (h)<br>14.6 in (w) x 21.3 in (d) x 10.2 in (h)<br><br>Weight: 16 kg (35.3 lb)   |
| Laser Test Specifications      | Class 1 Laser Product according to IEC 60825-1:2014 and EN 60825-1:2014<br>Complies with FDA performance standards for laser products 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007. |
| External Fusing                | Main fuse (appliance inlet): 4 A, 250 V; IEC 60127 rated as T4AL 250 V   |

Table 1-1: Specifications

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**Table 1. Preparation of High Controls**

| NUMBER OF UNIQUE CONTROL CONCENTRATES | VOLUME OF EACH CONTROL CONCENTRATE TO ADD TO MAKE HIGH CONTROL (μL) | VOLUME OF SAMPLE DILUENT TO ADD TO MAKE HIGH CONTROL (μL) |
|---------------------------------------|---|---|
| 1                                     | 120   | 3880  |
| 2                                     | 120   | 3760  |
| 3                                     | 120   | 3640  |
| 4                                     | 120   | 3520  |
| 5                                     | 120   | 3400  |
| 6                                     | 120   | 3280  |
| 7                                     | 120   | 3160  |
| 8                                     | 120   | 3040  |

**Table 2. Preparation of Low Controls**

| VOLUME OF HIGH CONTROL OR SAMPLE DILUENT TO ADD TO MAKE LOW CONTROL (μL) |      |
|--|------|
| High Control   | 80   |
| SDXX   | 3920 |

(XX=Diluent type\*)

**Table 3a. Preparation of BAFF/CCL8 High Controls**

| HIGH CONTROL DILUTION EXCEPTIONS | VOLUME OF CONTROL CONCENTRATE TO ADD TO MAKE HIGH CONTROL (μL) | VOLUME OF SAMPLE DILUENT TO ADD TO MAKE HIGH CONTROL (μL) |
|----------------------------------|--|---|
| BAFF                             | 30   | 3970  |
| CCL8                             | 40   | 3960  |

**Table 3b. Preparation of ITAC Low Controls**

| VOLUME OF ITAC HIGH CONTROL OR SAMPLE DILUENT TO ADD TO MAKE LOW CONTROL (μL) |      |
|---|------|
| High Control  | 40   |
| SD16  | 3960 |

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
**Bioanalytical qualification of clinical biomarker assays in plasma using a novel multi-analyte Simple Plex™ platform** *BioAnalytical Sciences, Oncology Biomarker Development Genentech Inc, South San Francisco Bioanalysis (2016) 8(23), 2415–2428*

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### Sample Preparation


Serum and plasma samples require a minimum 2-fold dilution with sample diluent. Samples above the ULOQ require further dilution.

### Sample Preparation

Serum and plasma samples require a minimum 10-fold dilution with sample diluent. Samples above the ULOQ require further dilution.

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




The screenshot shows the Simple Plex Explorer v.3.5.0.10 interface. At the top, there's a menu bar and a table with columns: Encrypted, Filename, Run Date, User, Specimen Type(s), Species Type, Lot Number, and Cartridge. Below this is a section for Kit Id: 79632 - Cartridge id: 79562, with fields for Analytes (OPN) and Comments. The main data table has columns for sample ID, channel, concentration, and other parameters. Some values are highlighted in pink (e.g., 56.13, 20501, 208.18, 94580) and some have asterisks (e.g., 329637\*, 166401). A dialog box is open for RFU %CV, showing values like 1.44 and 2.01. A settings menu is also open, showing options like 'Set Result Warning Highlighting To' and 'Use Filter Setting'. Red annotations explain that pink highlights indicate RFU CV above a threshold and black asterisks indicate concentrations outside the linear range.

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





**Load an existing cartridge configuration**

Please sign-in to load or save your selections. When you click the link, you'll be taken to the sign-in page, then brought back here to pick-up where you left-off.

Enter configuration code or 6-digit part number

**Make your selections below to configure your assay. Learn more about Simple Plex Assay Configurations**

**Species**

**Sample Type**

**Analytes**

Select from the list below to configure panels of up to 4 analytes.  
 Our team is here to help! Please reach out to your local Bio-Techne immunoassay specialist if you are having difficulty configuring your assay or are interested in analytes not seen in the menu below.

**Panel Builder Summary**

**Species:** Human  
**Sample Type:** Cell Supernatant


**Panel 1 Analytes:**  
 Diluent: SD13  
 Dilution Factor: 1:2  
 1. IL-6 2nd gen

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### Simple Plex Mouse/Rat TGF-beta 1 Cartridge

Catalog # SPCKB-MP-001191

Simple Plex Cartridge Kit containing TGF-b1 for use with Mouse and Rat Plasma/Serum. Includes wash buffer and sample diluent. For use with the Ella automated immunoassay system.

| Catalog #       | Availability | Size / Price                                    | Qty |
|-----------------|--------------|---|-----|
| SPCKA-MP-001191 | In Stock     | 1 Kit (1 Analyte x 16 Samples) / Please Inquire |     |
| SPCKC-MP-001191 | In Stock     | 1 Kit (1 Analyte x 32 Samples) / Please Inquire |     |
| SPCKB-MP-001191 | In Stock     | 1 Kit (1 Analyte x 72 Samples) / Please Inquire |     |

proteinsimple

[Datasheet / SDS](#)

**Order Information**

**Bio-Techne China Co. Ltd.**  
tel: 400-821-3475 (技术热线)  
info.cn@bio-techne.com

**Shipping Information**  
Please contact us for shipping charges and delivery times

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[specification sheet](#)


### Mouse/Rat TGF-β1, 2nd Generation

Simple Plex assay for the detection of mouse and rat Transforming Growth Factor beta 1 (TGF-β1) in serum and platelet-poor plasma (EDTA/Heparin)

For research use only. Not for use in diagnostic procedures.

#### Calibration Curve

The factory generated calibration curve shown below was compiled by averaging 5 replicates of each calibrator from multiple runs. The 5PL curve fit shows calibrator concentration as a function of signal intensity (relative fluorescent units, RFU).

#### Endogenous Levels

Endogenous levels were calculated from multiple samples.

| MOUSE SAMPLES           | MEAN (pg/mL) | RANGE (pg/mL) | STD DEV (pg/mL) |
|-------------------------|--------------|---------------|-----------------|
| Serum (n=5)             | 85,018       | 65,642-99,721 | 12,298          |
| PP EDTA Plasma (n=5)    | 10,130       | 5907-12,992   | 2645            |
| PP Heparin plasma (n=5) | 25,415       | 16,344-32,101 | 7593            |

| RAT SAMPLES             | MEAN (pg/mL) | RANGE (pg/mL) | STD DEV (pg/mL) |
|-------------------------|--------------|---------------|-----------------|
| Serum (n=5)             | 55,655       | 1491-85,661   | 32,823          |
| PP Heparin plasma (n=5) | 14,786       | 9394-18,669   | 3649            |

#### Precision

**Intra-assay Precision:** Each control was tested 16 times in one assay.

**Inter-assay Precision:** Replicates of each control were

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